



P-Sperm™ Video & Data Management Software for the SQA-Vp



Video & Data Management Software for use with the SQA-Vp

P-Sperm™ is a software package that is included with every SQA-Vp and greatly enhances the data management, storage and video capabilities of the system.

Features

- Import test data from the SQA-Vp
- Automatic AI dosing: Calculate doses by motile, total or morphologically normal sperm
- Collect and sort data, run graphs and view the test results by individual boar
- View live semen sample images on the PC screen
- Save test results with video clips or still images attached
- Password/security
- Optional language settings

System Requirements

- PC: 1 GHz processor, Pentium 3, RAM: 256 MB
- Monitor: 15" color
- AGP-video display card with at least 16 MB of RAM memory
- Video color: At least 16 bit (65,535)
- CD ROM drive
- 10 GB free hard disk space for image capturing (approx. 10,000 clips)
- Video resolution: Minimum 640 x 480
- Operating system compatibility: Windows XP, Windows Vista
- Ports: One serial; two USB ports
- MS Excel/Word