

# QwikLink™ LIS

Link to Laboratory Information Systems

*Transfer test and patient data from the SQA-V Gold to the Laboratory Information System*



The **QwikLink™ LIS** is a software program written specifically for the SQA-V *Gold* and V-Sperm *Gold* software to provide an interface to the Laboratory Information System. QwikLink™ LIS can be customized to transfer only the information selected by the user.

QwikLink™ LIS performs the following functions:

- Starts when the V-Sperm software is activated
- Runs as a background program while the SQA-V and V-Sperm are running
- Retrieves data from the V-Sperm database and sends to an external host

Each QwikLink™ LIS kit contains:

- One QwikLink™ LIS installation disk
- Two hours of tier two technical service support
- Full instructions for customized data transfer
- Complete list of database field descriptions
- Data flow overview
- RS232 converter and installation software
- Catalog #6400